




IONIC LIQUID GC COLUMNS

Extending the dimensions in Gas Chromatography

Anders Fridström
Technical Marketing Manager, Advanced Analytical



MERCK

Agenda

- Overview of Ionic Liquids
- New Ionic Liquid Column Developments
 - Watercol™ series
 - SLB® IL (i-series)
- Summary/Conclusions

Ionic Liquids

Properties as GC Columns

- A class of ionic solvents with low melting points
- Remain liquid over a wide temperature range (Room Temperature \rightarrow 350 °C)
- Unique combination of cations and anions that can provide different **selectivities** when used as stationary phases in GC
- Broadest range of solvation interactions of any known solvent



Unique Selectivity

The Importance of Selectivity

- A column's **selectivity** has the greatest influence on **resolution**




Ionic Liquids

- Extensive evaluations of ionic liquid GC columns
 - Main strength is **unique selectivity**
 - Often resulting in
 - Increased **resolution**
 - Shorter run times

Geminal Dicationic Ionic Liquid Stationary Phase


SLB-IL100

1,9-di(3-vinyl-imidazolium) nonane bis(trifluoromethyl) sulfonyl imidate



GC Column Polarity Scale


Visual Representation



Unique Selectivity


Fuel Alcohols; 90 °C Isothermal

SLB-IL60
30 m x 0.25 mm I.D., 0.20 µm




1. 2-Methyl-1-Butanol (Active Amyl Alcohol)
2. 3-Methyl-1-Butanol (Isoamyl Alcohol)

PEG
30 m x 0.25 mm I.D., 0.25 µm



Polarity vs. Selectivity

Principle Component Analysis (PCA) of



Solvation Parameter Model (SPM) Data

Only the ionic liquids are capable of simultaneously providing:

- ✓ Intense **H-acceptor interactions** (*a* constant)
- ✓ Intense **H-donor interactions** (*b* constant)
- ✓ **Dipolar interactions** (*s* constant)
- ✓ **π - π interactions** (*e* constant)
- ✓ Limited **dispersive interactions** (*l* constant)


95 probes / 52 columns

- Line 1: Poly(trifluoropropyl)siloxane
- Line 2: Poly(phenyl)siloxane
- Line 3: Poly(cyanopropyl)siloxane
- Cluster I: Polyethylene glycol
- Cluster II: Ionic liquid

Ref: Supelco Reporter 33.1 (March 2015), page 3-4.

Courtesy of Rosa Lebrón-Aguilar, CSIS, Madrid, Spain

Comparison of GC-MS TIC Bleed



All TICs are on the same Y-scale


SLB-IL100 GC-MS Bleed Profile

Abundance

Scan 6947 (24.821 min): 200804150006.D

MS spectrum of baseline at 230 °C

Different m/z bleed fragments compared to predominant m/z 207 and m/z 281 typically formed by polysiloxane polymers



column: SLB-IL100, 30 m x 0.25 mm I.D., 0.20 µm (28884-U)

oven: 60 ° C (1 min.), 8 ° C/min. to 230 ° C (5 min.)

inj.: 250 ° C

MSD interface: 220 ° C

scan range: m/z = 35-500

carrier gas: helium, 1.5 mL/min. constant

injection: 1 µL, splitless (1.0 min.)

liner: 4 mm I.D., single taper

Watercol™ Series Capillary GC Columns


Convenient Measurement of Water



Watercol Series Capillary GC Columns

Overview


- Contain innovative ionic liquid stationary phases
- Produce a sharp peak shape for water




Three chemistries, each with a different selectivity.

Watercol™ Series


Watercol™1460



Watercol™1900



Watercol™1910



Watercol Series Capillary GC Columns

Key Benefits of a Sharp Water Peak Shape

- **Measurement of water** (qualitative and quantitative)
- **Linear response** over a very wide range (0.01% water to 100% water)
- **Great sensitivity**
 - 100 ppm using a thermal conductivity detector (**TCD**)
 - Anticipate much better sensitivity with
 - A mass spectrometer (**MS**) operated in the SIM mode
 - Vacuum ultraviolet (**VUV**) absorption spectroscopy
 - A barrier discharge ionization detector (**BID**)
- **Reproducibility** (virtually no change in column performance over time)


NOTE: The detector being used must be able to detect water.

Watercol Series Capillary GC Columns


Selectivity Options; 96 °C, using a GC/TCD

1. n-tridecane
2. 2-octanone
3. 1-heptanol
4. 1-octanol
5. naphthalene


Watercol™1460



Watercol™1900



Watercol™ 1910




Comparison of water selectivity and peak shape

GC/TCD, Temp. 80 °C (isothermal)


Watercol™ 1910

30-m x 0.25 mm i.d x 0.20 µm df.




Wax-type column

30-m x 0.25 mm i.d x 0.25 µm df.




Q-Plot

30-m x 0.32 mm i.d




Water Calibration Curve (0.01-0.5%) on Watercol™1910




Watercol Series Capillary GC Columns

Measure Water in a 60-Component Solvent Mix (Chemical Application)


- Water peak is baseline resolved from a complex matrix that consists of the most common solvents used in chemical/API production
- Good peak symmetry, suitable for water quantification



GC/FID




GC/TCD



Watercol Series Capillary GC Columns

Measure Water in Gasoline (Petroleum Application)

- Water peak is baseline resolved from other gasoline components
- Analysis of BTEX cmpds, and oxygenates such as Ethanol are also possible. Multiple target analytes can be identified in a single run



Watercol Series Capillary GC Columns

Comparison to Karl Fischer

- Karl Fischer titration (mostly used)

- 500 ppm is typical detection limit
- **Watercol + TCD = 100 ppm**

- Karl Fischer coulometric


- 5 ppm is typical detection limit
- Watercol + other detectors =

Watercol columns have several advantages over Karl Fischer

- Reduced **costs**
 - GC requires **less hands-on time**
 - Lower volume of **chemical waste** generated (which must be disposed of)
 - Perform two analysis (water and volatiles/semivolatiles) with a **single method**
- Reduced **amount of sample** needed (<0.5 mL for liquids and <0.2 g for solids)
- **Eliminate worker exposure** to (potentially) **harmful chemicals**
- **No unwanted side reactions**
- **No solubility issues**

Reproducibility and Detection Limits VUV

A



B

Parameter	Setting	Unit
Average Area	0.00187	AU
% RSD	4.68	-
Conc. Background	100	ppm
MDL	13	ppm
Average Height	0.035	AU
RMS Noise	.00029	AU
S/N	120	-
LOD	2.5	ppm

Shimadzu BID Detector




Fig. 1: Chromatogram for water determination (25 ppm) in LPG. Quantification of Limit(S/N=10) and Detection of Limit(S/N=3.3) can be down to 0.66 ppm and 0.22 ppm respectively.

Watercol 1910

Experimental Run Conditions

Instrument:	Shimadzu QP-2010 Plus GC
Column:	Watercol 1910, 30 m x 0.25 mm I.D., x 0.2 µm df
Oven:	55°C (2min), 4°C/min to 180°C (1 min)
Inj. Temp:	250°C
Carrier Gas:	Helium, (40 cm/sec)
Detector:	BID, 200°C
Liner:	4 mm I.D. split liner with wool
Injection:	0.5 µL, 200:1 split
Sample:	SPME Alcohol Mix in H ₂ O

Analysis of water in Liquid Petroleum Products





Fig. 1: Chromatogram for water determination (25 ppm) in LPG. Quantification of Limit($S/N=10$) and Detection of Limit($S/N=3.3$) can be down to 0.66 ppm and 0.22 ppm respectively.


Water as Solvents

Alcohols and other VOC analysis in water; 35 °C, 4 °C/min to 125 °C (2 min)



Watercol Series Capillary GC Columns

Measure Aroma Fusel Alcohols in Tequila (Food and Beverage Application)



Watercol Series Capillary GC Columns

Summary

- Measurement of water with linearity, great sensitivity, and reproducibility
- Less hands-on time compared to other technologies used to measure water
- Less volume of chemical waste
- Possible to obtain results for water + volatiles/semivolatiles in same analysis
- Amendable for water as injection solvent

GC columns made with ionic liquid stationary phases are exclusive to Supelco.

SLB[®]-IL (i-series) Capillary GC Columns

Polar Selectivity and Inertness



SLB-IL (i-series) Capillary GC Columns

Overview

- More inert versions of popular ionic liquid chemistries
- Provide selectivity **and** inertness for **polar analytes**
- Polar selectivity and inertness allows
 - **Better accuracy** (identification and amount)
 - **Better sensitivity**

SLB-IL (i-series) Capillary GC Columns


Selectivity Options

- A range of i-series columns were developed and classified as
 - Polar (SLB-IL60i)
 - Highly polar (SLB-IL76i)
 - Extremely polar (SLB-IL111i)
 - **Great choices for 2D GC**

sigma-aldrich.com/il-gc-inert

Ionic Liquids GC columns are great columns for GC×GC

Increase orthogonality when paired with a non-polar stationary phase (e.g. SLB-5ms)




Courtesy of Prof. L. Mondello

SLB-IL 60i

Results-Xylenols Mix 130°C Isothermal


Peak	Compound	t _R
1	2,6-xlenol	4.825
2	2,4-xlenol	7.019
3	2,5-Xylenol	7.236
4	2,3-Xylenol	8.447
5	3,5-Xylenol	10.317
6	3,4-Xylenol	11.769



Industrial Solvents

Figure 2. Industrial Solvents on SLB-IL60i

column: SLB-IL60i, 30 m x 0.25 mm I.D., 0.20 µm (29832-U)
 oven: 40 °C (4 min), 8 °C/min to 200 °C (5 min)
 inj. temp.: 250 °C
 detector: FID, 250 °C
 carrier gas: helium, 30 cm/sec
 injection: 1 µL, 100:1 split
 liner: 4 mm I.D., split type, cup design
 sample: industrial solvents, each at 0.2 % (w/v) in pentane



- | | |
|-------------------------------|--|
| 1. Hexane | 29. 1,2-Dichloroethane |
| 2. 1,1-Dichloroethylene | 30. <i>n</i> -Butyl acetate |
| 3. Methyl formate | 31. 2-Hexanone |
| 4. Acetone | 32. Isobutanol |
| 5. Ethyl formate | 33. Nitropropane |
| 6. Methyl acetate | 34. Isoamyl acetate |
| 7. trans-1,2-Dichloroethylene | 35. Ethylbenzene |
| 8. Tetrahydrofuran | 36. Mesityl oxide |
| 9. Carbon tetrachloride | 37. <i>p</i> -Xylene |
| 10. 1,1-Dichloroethane | 38. <i>m</i> -Xylene |
| 11. Ethyl acetate | 39. 5-Methyl-2-hexanone |
| 12. Methanol | 40. <i>n</i> -Butanol |
| 13. Isopropyl acetate | 41. <i>n</i> -Amyl acetate |
| 14. 2-Butanone | 42. <i>o</i> -Xylene |
| 15. 2-Propanol | 43. Isoamyl alcohol |
| 16. Methylene chloride | 44. Chlorobenzene |
| 17. Ethanol | 45. Styrene |
| 18. Benzene | 46. 1,1,1,2-Tetrachloroethane |
| 19. <i>n</i> -Propyl acetate | 47. Dimethylformamide |
| 20. Trichloroethylene | 48. Diacetone alcohol |
| 21. 4-Methyl-2-pentanone | 49. Cyclohexanol |
| 22. Isobutyl acetate | 50. 2-Butoxyethanol (Butyl cellosolve) |
| 23. Tetrachloroethene | 51. 1,4-Dichlorobenzene |
| 24. Chloroform | 52. 1,1,2,2-Tetrachloroethane |
| 25. sec-Butanol | 53. 2-Methylphenol |
| 26. Toluene | 54. 3-Methylphenol |
| 27. <i>n</i> -Propanol | 55. 4-Methylphenol |
| 28. 1,4-Dioxane | |

Industrial Solvents

Selectivity Difference 111i and 60i

column: SLB-IL111i, 30 m x 0.25 mm I.D., 0.20 µm (29883-U)

oven: 40 °C (8 min), 8 °C/min to 200 °C (1 min)

All other conditions and peak IDs are the same as Figure 2.




Figure 2. Industrial Solvents on SLB-IL60i

column: SLB-IL60i, 30 m x 0.25 mm I.D., 0.20 µm (29832-U)

oven: 40 °C (4 min), 8 °C/min to 200 °C (5 min)

inj. temp.: 250 °C


detector: FID, 250 °C

carrier gas: helium, 30 cm/sec

injection: 1 µL, 100:1 split

liner: 4 mm I.D., split type, cup design


sample: industrial solvents, each at 0.2 % (w/v) in pentane



- | | |
|---------------------------------------|--|
| 1. Hexane | 29. 1,2-Dichloroethane |
| 2. 1,1-Dichloroethylene | 30. <i>n</i> -Butyl acetate |
| 3. Methyl formate | 31. 2-Hexanone |
| 4. Acetone | 32. Isobutanol |
| 5. Ethyl formate | 33. Nitropropane |
| 6. Methyl acetate | 34. Isoamyl acetate |
| 7. <i>trans</i> -1,2-Dichloroethylene | 35. Ethylbenzene |
| 8. Tetrahydrofuran | 36. Mesityl oxide |
| 9. Carbon tetrachloride | 37. <i>p</i> -Xylene |
| 10. 1,1-Dichloroethane | 38. <i>m</i> -Xylene |
| 11. Ethyl acetate | 39. 5-Methyl-2-hexanone |
| 12. Methanol | 40. <i>n</i> -Butanol |
| 13. Isopropyl acetate | 41. <i>n</i> -Amyl acetate |
| 14. 2-Butanone | 42. <i>o</i> -Xylene |
| 15. 2-Propanol | 43. Isoamyl alcohol |
| 16. Methylene chloride | 44. Chlorobenzene |
| 17. Ethanol | 45. Styrene |
| 18. Benzene | 46. 1,1,1,2-Tetrachloroethane |
| 19. <i>n</i> -Propyl acetate | 47. Dimethylformamide |
| 20. Trichloroethylene | 48. Diacetone alcohol |
| 21. 4-Methyl-2-pentanone | 49. Cyclohexanol |
| 22. Isobutyl acetate | 50. 2-Butoxyethanol (Butyl cellosolve) |
| 23. Tetrachloroethene | 51. 1,4-Dichlorobenzene |
| 24. Chloroform | 52. 1,1,2,2-Tetrachloroethane |
| 25. <i>sec</i> -Butanol | 53. 2-Methylphenol |
| 26. Toluene | 54. 3-Methylphenol |
| 27. <i>n</i> -Propanol | 55. 4-Methylphenol |
| 28. 1,4-Dioxane | |

SLB-IL (i-series) Capillary GC Columns

Active Amyl Alcohol and Isoamyl Alcohol with n-Alkanes; 90 °C



Applications

- Critical isobaric PAHs
- TCDD
- Allergens
- Industrial Solvents
- FAMES
- BTEX
- and Much more..

sigma-aldrich.com/GC

CONCLUSIONS

- **Something totally new and completely different in the world of GC phases**
- **Have the opportunity to impact current GC and GC-MS practices along several paths**

Watercol™ Series

Watercol 1460

Watercol 1900

Watercol 1910

Three unique different selectivity ionic liquid phases for the GC analysis of H₂O

“I” Series column provide unique selectivity with improved inertness

sigma-aldrich.com/gc

Web Site (*sigma-aldrich.com/gc*)



Analytical / Chromatography Products

+ Air Monitoring

+ Analytical Reagents

+ Analytical Standards

Applications Search

Bioanalysis

+ Chiral Chromatography

Custom Products

+ Food and Beverage Analysis

- Gas Chromatography

Column Selection

Gas Chromatography Products

+ GC Learning Center

+ Capillary GC Columns

Packed GC Columns

GC Column Test Mixes

GC Accessories

Gas Chromatography

Supelco's full-line of capillary gas chromatography (GC and GC-MS) columns includes ionic liquid, MS-grade, Fast GC, GCxGC, chiral, and PLOT columns, conveniently grouped by industry (including environmental, petroleum, food & beverage, and pharmaceutical), by method reference, by phase polarity and also by name. We offer the most complete line of packed GC columns and components. To maximize performance, check out our comprehensive line of GC accessories (including molded GC septa, deactivated inlet liners, syringes, and vials), standards, solvents, and gas purification/management products (including purifiers, gas generators, tubing, fittings, and pressure regulators).

Products:

- Gas Chromatography Products
- Capillary GC Columns and Guard Columns
- Ionic Liquid GC Columns
- Packed GC Columns and Components
- GC Column Test Mixes
- GC Accessories
- Purge and Trap
- Analytical Standards
- GC Solvents
- GC Derivatization Reagents
- Gas Purification/Management
- Adsorbents

Tool Box:

- Interactive Supelco Catalog
- Custom Capillary GC Column Request **NEW**
- Custom Packed GC Column Request **NEW**

Special Offers:

- GC Special Offers
- Sigma-Aldrich Monthly Savings

Analytical / Chromatography

Learning Center:

- Column Selection
- Chromatograms
- Presentations
- Brochures
- Literature

Related Topics:

- Fast GC
- GCxGC
- Environmental Testing
- Dioxin Sample Prep
- Biofuels Manufacturing & Testing
- Food & Beverage Analysis
- Fatty Acid & FAME Analysis
- Solid Phase Microextraction (SPME)

