

Sample Preparation Techniques and Method Development

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Sep-Pak®
Sample Extraction Products



Sample Preparation

DISQUE™
Dispersive Sample Preparation

Waters
ASIS®
SAMPLE EXTRACTION PRODUCTS

Ostro™
SAMPLE PREPARATION PRODUCTS

Sirocco™
Protein Precipitation Plate

PoraPak™
Rxn
Post Synthesis Cleanup

Contact:

evelien_dejaegere@waters.com
marleen_van_wingerden@waters.com

Our Approach to Method Development

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- No “one size fits all”
- Different segments of development process
 - Scientific and business drivers may be different
 - Drivers may be the same but with varying degrees of risk tolerance
- Use of scientifically appropriate criteria for final method choice

Overview

- Food en Environmental



- Bioanalysis



- Synthesis



Dispersive Sample Preparation

"QuEChERS" Method

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™


- "Quick, Easy, Cheap, Effective, Rugged, Safe"
- Popular approach to performing sample cleanup
 - Low material costs
 - Method gaining world-wide acceptance
- High throughput sample preparation and screening analysis
- Combination of salting-out LLE and matrix dispersion extraction

DisQUE™
Dispersive Sample Preparation

Dispersive Sample Preparation *Benefits*

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- Low cost
- Minimum specialist skill
- Reduced sample sizes
- Low development and troubleshooting



DISQUE™
Dispersive Sample Preparation

Typical Laboratory Workflow

- Most “QuEChERS” and dispersive SPE applications require the user to source the individual components including vials, sorbents and buffers.
 - DisQuE Dispersive Sample Preparation Kits eliminate this lab outsourcing
 - Improves consistency in extraction product
 - Improves lab efficiency and workflow
 - Single Vendor source improves quality

- DisQuE offers:
 - Cleaner devices with lower extractables
 - Method and product support to enable customer success

Extraction Tubes (Tube 1)

<i>AOAC Configuration</i>	<i>CEN Configuration</i>
1.5 g Sodium Acetate 6 g Magnesium Sulphate	4 g Magnesium Sulphate 1 g Sodium Chloride 1 g Trisodium Citrate 0.5 g Disodium Citrate

Clean Up Tubes (Tube 2 - 2 mL Option)

<i>AOAC Configuration</i>		<i>CEN Configuration</i>	
150mg Magnesium Sulphate 50mg PSA	150mg Magnesium Sulphate 50mg PSA 50mg C ₁₈	150mg Magnesium Sulphate 25mg of PSA	150mg Magnesium Sulphate 25mg of PSA 25mg C ₁₈

Clean Up Tubes (Tube 2 - 15 mL Option)

900 mg Magnesium Sulphate 150 mg of PSA	900 mg Magnesium Sulphate 150 mg of PSA 150 mg C ₁₈
--	--


Extraction Tube 1

- Water is IMPORTANT!!
 - Removes unwanted matrix
 - Acetonitrile extraction more effective

Acetonitrile Layer (**ANALYTES ARE HERE**)

Water Layer (**POLAR INTERFERENCES**)

Solids Layer (**BUFFER SALTS AND SAMPLE**)




Cleanup Tube 2

- Provides additional cleanup
- Sorbent choices
 - PSA
 - Graphitized Carbon Black (GCB)
 - C18

Acetonitrile Layer (**ANALYTES ARE HERE**)

Sorbent (**INTERFERENCES**)



DisQuE™ Performance Highlights from Application Note

Waters
THE SCIENCE OF WHAT'S POSSIBLE™

- DisQuE method tested against 402 pesticide residues
 - Application Note: 720002628en
- Approved AOAC Method
 - Ref. Lehotay, *JAOAC Int.* **90**(2) 185-520 (2007).
- Excellent recovery
- Consistent results from matrix to matrix

[APPLICATION NOTE]

Waters
THE SCIENCE OF WHAT'S POSSIBLE™

A RAPID METHOD FOR THE SCREENING AND CONFIRMATION OF OVER 400 PESTICIDE RESIDUES IN FOOD

James Morphet and Peter Hancock, Waters Corporation, Manchester, UK

AIM

To utilize the power of UltraPerformance Liquid Chromatography (UPLC®) combined with fast MS acquisition rates, to give a rapid method for the screening of 402 pesticide residues in a single 10 minute run. A second injection, for confirmatory purposes, will meet SANCO Analytical Quality Control procedures for pesticide residue analysis (SANCO/2007/31311).

Advances in chromatographic separation and detection technologies have enabled analysts to increase the number of analytes determined in a single run. Tandem quadrupole mass spectrometry offers a highly specific and selective detection technique that has become the technique of choice within the laboratory.²

The following method describes a solution for the rapid analysis of pesticides in mango, avocado, and fruit based baby food extracts that is able to exceed current worldwide legislation.

INTRODUCTION

Pesticides are widely used in the production of foodstuffs to meet consumer demand for plentiful food at reasonable prices, all year round. However, continued growth in the use of pesticides, poor agricultural practices, and illegal use can pose significant risks to human health through the presence of pesticide and metabolite residues in food products. Most countries have strict regulations governing pesticides. Legislation imposes Maximum Residue Limits² (MRLs) for pesticide residues in food products requiring analytical techniques that are sensitive, selective, and robust.

Multi-residue pesticide analysis is challenging due to the low levels present, the wide variety of pesticides, and the very different chemical classes they represent. As there are currently well over 1,000 pesticides in use, laboratories are under increasing pressure to broaden the range of pesticides determined in a single analysis over a shortened run time.

The need to meet mandated detection limits, develop generic sample preparation techniques for complex matrices, and the desire to increase sample throughput are the main challenges facing food safety testing laboratories today. The use of a single multi-residue method per instrument can dramatically improve return on investment by removing the need to change method parameters. This is often the case when analysing a wide variety of commodities with differing lists of legislated pesticides.

EXPERIMENTAL

A QuEChERS extraction was utilized for this multi-residue method by homogenizing a food sample, adding organic solvent, salts, and sorbent. Once mixed, the pesticide residues were partitioned into the organic solvent, which was then subjected to further clean up. The supernatant was collected, diluted, and injected onto the LC/MS/MS system as described below.

Extraction Procedure⁴:

1. Add 15 g homogenized sample to a 50 mL extraction tube containing 1.5 g sodium acetate and 6 g magnesium sulfate. Add 15 mL 1% acetic acid in acetonitrile.
2. Add any pre-extraction internal standards.
3. Shake vigorously for one minute and centrifuge > 1500 rcf for one minute.
4. Transfer 1 mL of the acetonitrile extract into the 2 mL centrifuge tube containing 50 mg PSA and 150 mg of magnesium sulphate.
5. Shake for 30 seconds and centrifuge > 1500 rcf for one minute.
6. Transfer 100 µL of final extract into an autosampler vial. Add any post-extraction internal standards. Dilute with 900 µL water.

- Quechers




- Bioanalysis



- Synthesis



Different Methods for Different Purposes: Decision Making Process



Highest Selectivity



Direct Inject



Simple



High throughput



PL Removal



Highest Sensitivity



"Clean" Extracts




Typical Sample Preparation Selection Process

1. The simplest method which meets the assay needs is usually chosen.
2. For very challenging assays (low detection limits, closely related endogenous constituents, inhalation products, peptides, etc), SPE is often the first choice.
3. Exact technique chosen will depend on outcome of study and how much risk can be tolerated.


Method 1: Clopidogrel and Ticlopidine (IS) in Plasma

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Clopidogrel
MW 321.8
pKa = 4.5



Ticlopidine (IS)
MW 263.8
pKa = 7.6



Assay Use:

- Routine analysis of patient samples, GLP or clinical lab

Assay Requirements:





- LLOQ 10-25 pg/mL
- Simple method
- Must transfer across lab with varying levels of expertise
- Concerns about build up of phospholipids on LC columns and in MS source

Clopidogrel and IS: Chromatographic Screening Results





Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

ACQUITY UPLC® BEH C₁₈ Column

Clopidogrel




IS




Sample Preparation Comparison

% Analyte Recovery



% Matrix Effects




1. Adequate recovery?
2. Very simple method?
3. Phospholipid removal?

Benefits of Phospholipid Removal Plates for this Application

- Very simple protocol
 - Easily transferred from lab to lab
- Removes vast majority of phospholipids
 - Improved instrument uptime
 - More robust methods
- High recovery
 - Helps meet detection limits

Representative Validation Results: Ostro™ Sample Preparation Plate

Compound name: Clopidogrel
 Correlation coefficient: $r = 0.999471$, $r^2 = 0.998943$
 Calibration curve: $0.049871 \cdot x + 0.000843056$
 Response type: Internal Std (Ref 2), Area * (IS Conc./IS Area)
 Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



	Standard conc. ng/mL	Clopidogrel Area	IS Area	Response	Calc. conc. ng/mL	%Dev
Standard	0.01	157.3	122606.4	0.001	0.009	-11.8
Standard	0.05	426.5	118007.4	0.004	0.056	11.1
Standard	0.1	560.7	104434.5	0.005	0.091	-9.2
Standard	0.5	3013.1	112846.8	0.027	0.518	3.7
Standard	1	6155.2	114711.6	0.054	1.059	5.9
Standard	5	28044.7	114521.0	0.245	4.894	-2.1
Standard	10	60742.4	115491.1	0.526	10.529	5.3
Standard	50	274997.2	115583.9	2.379	47.690	-4.6
Standard	100	590079.4	116193.5	5.078	101.814	1.8
QC	0.075	432.0	109556.5	0.004	0.062	-17.1
QC	0.75	4382.0	113846.2	0.038	0.755	0.7
QC	7.5	42802.8	113869.8	0.376	7.520	0.3
QC	75	447395.8	110303.7	4.056	81.314	8.4
QC	0.075	355.6	112084.7	0.003	0.047	-37.7
QC	0.75	4944.7	121039.7	0.041	0.802	7
QC	7.5	39973.1	110003.2	0.363	7.270	-3.1
QC	75	403244.8	108186.3	3.727	74.722	-0.4
QC	0.075	534.3	123902.3	0.004	0.070	-7.2
QC	0.75	5479.0	126105.7	0.043	0.854	13.9
QC	7.5	46733.5	118830.3	0.393	7.869	4.9
QC	75	317429.4	93421.1	3.398	68.116	-9.2

Assay Requirements:

- LLOQ 10-25 pg/mL
- Simple method
- Must transfer across lab with varying levels of expertise
- Lab is very concerned about build up of phospholipids on LC columns and in MS source

Sample Preparation Options

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

SIMPLY CLEAR

THE FASTER WAY TO CLEANER

SENSITIVITY IN ITS PUREST FORM

The image displays three Waters sample preparation products against a dark blue background. On the left is a Sirocco Protein Precipitation Plate, a 96-well plate with a white label that includes the Waters logo, 'Sirocco™ Protein Precipitation Plate', 'PIN: 196002448', and 'Lot No. 906801'. In the center is an Ostro 96-Well Plate, a clear 96-well plate with a white label that includes the Ostro logo, 'Waters THE SCIENCE OF WHAT'S POSSIBLE.™', 'Ostro™ 96-Well Plate', 'Mfg. Lot No.: 123456789B', 'Sorbent Batch No.: 101', 'Part No.: 196005518', and 'Made in Ireland'. On the right is an Oasis HLB μElution Plate, a 96-well plate with a white label that includes the Oasis logo, 'Oasis® HLB μElution Plate', 'Part Number 1960016288A', and 'Lot Number 057A35215A'. Below the plates are the Sirocco Protein Precipitation Plate logo, the Ostro SAMPLE PREPARATION PRODUCTS logo, and the Oasis SAMPLE EXTRACTION PRODUCTS logo.

Sirocco™
Protein Precipitation Plate

Ostro™
SAMPLE PREPARATION PRODUCTS

OASIS®
SAMPLE EXTRACTION PRODUCTS

NEW Ostro™ Sample Preparation Plate

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Simple Pass-through Method

Basic Protocol

Ostro™
SAMPLE PREPARATION PRODUCTS



Add plasma

Add organic


Aspirate to mix

Apply vacuum

Analyze


Ostro™ Sample Preparation Plate: Analyte Recoveries

Best Recovery for Diverse Analytes



Ostro™ Phospholipid Removal: Reproducibility

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™





Reproducibility of Phospholipid Removal
in 6 Lots of Human Plasma

Method 2: Bioequivalence Study

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Ropinirole
MW 260.38
pKa 9.7



IS: Citalopram
MW 324.40
pKa 9.5

Assay Use:


- Routine analysis of patient samples, GLP or clinical lab

Assay Requirements:


- LLOQ 5 pg/mL
- Concerned with matrix effects
- Meeting FDA requirements
- Fast method development

Oasis® Solid-Phase Extraction (SPE)

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Oasis[®] 2x4 Method



SPE 96-Well Plate Format: Oasis® μ Elution Plates



Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- Novel plate technology enables 25 μ L SPE elution
 - Allows loading sample volumes from 25 to a maximum of 375 μ L
 - 50 to 750 μ L 1:1 diluted sample, 750 μ L is the well volume
 - Elution volume in as little as 25 μ L
- No evaporation means higher throughput and sensitivity
 - Sensitive and selective SPE for bioanalytical clinical samples
 - Increased sensitivity: up to 15x concentration factor (through format change)
- SPE without an evaporation step




Sample Enrichment: Up to a 15X Concentration Factor

0.5 ng/mL risperidone
375 μ L sample diluted with 375 μ L 4% H_3PO_4 loaded onto Oasis® MCX μ Elution plate



Improving Sensitivity When Sample Limited: Example

1 ng/mL Terfenadine in rat plasma, Oasis® MCX



Why Eliminate Evaporation and Reconstitution?

- Many general protocols include an evaporation and reconstitution step before injection
 - To obtain concentration necessary to reach desired limits of detection
- Certain small analytes such as ibuprofen and pseudoephedrine are susceptible to evaporative loss[†]
- Direct injection from the final eluate ensures maximum recovery for volatile analytes
- Direct injection saves analyst's time
- Direct injection eliminates an additional handling step
- Allows robots to be used to process more samples

[†]Naidong, W. *J Chromatog. B.* 796 (2003) 209-224


SPE Protocol: Oasis[®] MCX μ Elution 96-Well Plate

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Condition: 200 μ L MeOH
Equilibrate: 200 μ L H₂O
Load: Various volumes of plasma and urine,
diluted 1:1 with 4% H₃PO₄ in H₂O
Wash 1: 200 μ L 2% CHOOH in H₂O
Wash 2: 200 μ L MeOH
Elute: 2 x 25 μ L 5% NH₄OH in 60:40 ACN:MeOH
Dilute: 50 μ L H₂O
Inject: 5 μ L


Capacity: SPE Recovery for Polar and Non-Polar Analytes in Plasma Example

Oasis® MCX uElution Plate
200 ng/mL Imipramine and 200 ng/mL Atenolol




Capacity: SPE Recovery for Polar and Non-Polar Analytes in Urine Example

Oasis® MCX uElution Plate 200 ng/mL Imipramine and 200 ng/mL Atenolol



Matrix Effects for Polar and Non-Polar Analytes in Plasma

Oasis® MCX uElution Plate 200 ng/mL Imipramine and 200 ng/mL Atenolol





Loaded Plasma Volumes (Undiluted)

Method 2: Bioequivalence Study

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Ropinirole
MW 260.38
pKa 9.7



IS: Citalopram
MW 324.40
pKa 9.5

Assay Use:

- Routine analysis of patient samples, GLP or clinical lab

Assay Requirements:

- LLOQ 5 pg/mL
- Concerned with matrix effects
- Meeting FDA requirements
- Fast method development

SPE Protocol: Oasis[®] MCX μ Elution 96-Well Plate

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Condition: 200 μ L MeOH
Equilibrate: 200 μ L H₂O
Load: 600 μ L diluted Plasma sample (300 μ L plasma,
diluted 1:1 with 4% H₃PO₄ in H₂O)
Wash 1: 200 μ L 2% CHOOH in H₂O
Wash 2: 200 μ L MeOH
Elute: 25 μ L 5% NH₄OH in MeOH*

Inject: 8 μ L

LLOQ: 0.005 ng/mL

*During method development, it was determined that a 100% MeOH elution solvent was sufficient for good recovery

Ropinirole: Accuracy Data

FDA Acceptance Criteria: $\pm 15\%$, LLOQ $\pm 20\%$

Name	RT	Area	IS Area	Response	ng/mL	%Dev
blank plasma	1.77	32.8	70.0	0.468		
0.005 ng/mL	1.76	348.6	18914.4	0.018	0.005	-0.8
0.010 ng/mL	1.76	755.8	25514.2	0.03	0.01	0.3
0.020 ng/mL	1.76	1502.8	28613.7	0.053	0.02	1.9
0.100 ng/mL	1.76	6850.6	29148.2	0.235	0.103	3
1 ng/mL	1.76	60103.1	26299.2	2.285	1.031	3.1
5 ng/mL	1.76	267985.3	26320.0	10.182	4.604	-7.9
10 ng/mL	1.76	551196.6	24835.3	22.194	10.041	0.4

Deviations from -7.9% to 3.1%

Assay Requirements:



- Rapid method developed using Oasis 2x4 approach and Oasis μ Elution Plate



- LLOQ 5 pg/mL



- Method easily meets the FDA criteria for acceptability
 - Linearity across 3.5 orders of magnitude
 - Excellent sensitivity: 0.005 ng/mL (5 pg/mL)
 - Accuracy for all points on standard curve



- Matrix effects for analyte and internal standard <15%

SIMPLY CLEAR

THE FASTER WAY
TO CLEANER

SENSITIVITY IN ITS
PUREST FORM



Sirocco™
Protein Precipitation Plate



Ostro™
SAMPLE PREPARATION PRODUCTS



OASIS®
SAMPLE EXTRACTION PRODUCTS

- Quechers



- Bioanalysis



- Synthesis



Cleanup Challenges

- Convert from a high boiling point solvent before next reaction step
- Separate product (bases) from reactants and reaction solvents
- Remove TFA
- Convert from last solvent used in synthesis before Prep LC
- Remove water from fraction collected from Prep LC for faster dry down

Copolymer that exhibit the following properties:

- Hard material that does not develop increasing back pressure with flow. (Not compressing)
- Little swelling or shrinking across a range of solvents and pH extremes.
- Low hydraulic resistance enables flow by gravity.
- pH extreme tolerant without dissolution or hydrolysis, both limitations of silica-based sorbents.

This combination of physical and chemical properties makes Porapak Rxn an ideal material for synthesis cleanup.

Recommended Techniques



Reaction Products

	Technique	Recommended Product
Bases	Catch and Elute	PoraPak Rxn CX
Acids	Pass Through	PoraPak Rxn CX
Neutrals	Adsorption	PoraPak Rxn RP

Reaction Components

TFA removal	Catch and Elute	PoraPak Rxn CX
Convert from high boiling point solvent for bases	Catch and Elute	PoraPak Rxn CX

Pre-and-Post Preparative Column

Fractionate sample before Prep for bases	Catch and Elute	PoraPak Rxn CX
Remove water from post column fraction for bases	Catch and Elute	PoraPak Rxn CX

Sample Preparation Techniques and Method Development – Information?

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Product information on our website:

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

[Oasis – Sample Extraction Product](http://www.waters.com/Oasis)

<http://www.waters.com/Oasis>

[Ostro – Sample Preparation Product](http://www.waters.com/Ostro)

<http://www.waters.com/Ostro>

[Sirocco – Protein Precipitation Plate](http://www.waters.com/Sirocco)

<http://www.waters.com/Sirocco>

[Certified Sep-Pak – Critical Clean SPE](http://www.waters.com/waters/nav.htm?cid=10105548)

<http://www.waters.com/waters/nav.htm?cid=10105548>

[Sep-Pak – SPE](http://www.waters.com/SepPak)

<http://www.waters.com/SepPak>

[DisQue - Quechers](http://www.waters.com/DisQue)

<http://www.waters.com/DisQue>

[PoraPak – Synthesis Clean-Up](http://www.waters.com/PoraPak)

<http://www.waters.com/PoraPak>

[Envirogel GPC Clean-UP](http://www.waters.com/webassets/cms/support/docs/WAT036556.pdf)

<http://www.waters.com/webassets/cms/support/docs/WAT036556.pdf>